



Tempest Surgical Device Washer

Multiple mechanical actions to clean multiple device styles

Class II Device cleans, disinfects, and flushes cannulated surgical instruments such as:

- Suction tubes
- Long lumens with narrow diameters
- 24 ported devices in each lab tested cycle



Features:

- Efficient disinfection
- Enzymatic, agitated submersion
- Digital ultrasonic cleaning beyond OEM device requirements
- Controlled enzymatic seeding of interior lumen
- Multiple air injected interior lumen fluid flush
- High heat interior and exterior disinfecting rinse
- Final ozonated fluid spray
- Repeatable cleaning allows confident sterilization

For use in:

- Central Processing
- Ambulatory Surgical Centers
- OR Utility Rooms

Assure consistent, repeatable results that hand cleaning or ultrasonics alone cannot duplicate